



big squid IKAMAG® White

Magnetic stirrer in five motifs.

- Diameter: 160 mm
- Stirring quantity (H₂O) up to 1.5 l
- Recyclable materials
- Very good chemical resistance due to glass top and synthetic bottom made of Hytrel®
- Electronically controlled motor with infinitely variable speed from 0 - 1,500 rpm

Accessories: TRIKA® 25 Magnetic stirring bar, RSE Stirring bar remover, IKAFLON® 15 Magnetic stirring bar, IKAFLON® 25 Magnetic stirring bar, IKAFLON® 10 Magnetic stirring bar, IKAFLON® 20 Magnetic stirring bar, IKAFLON® 30 Magnetic stirring bar

Technical Data

Number of stirring positions	1
Stirring quantity max. per stirring position (H ₂ O) [l]	1.5
Stirring quantity max. (H ₂ O) [l]	1.5
Motor rating input [W]	2
Motor rating output [W]	1
Speed range [1/min]	0 - 1500
Stirring bar length max. [mm]	30
Set-up plate material	glass
Set-up plate dimensions [mm]	Ø 160
Dimensions (W x H x D) [mm]	180 x 48 x 180
Weight [kg]	1
Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80
Protection class according to DIN EN 60529	IP 54
Voltage [V]	100 - 240
Frequency [Hz]	50/60
Power input [W]	4
Ident. No.	3050009