

Data Sheet

Mixing / Magnetic Stirrers



topolino mobil

Extremely convenient magnetic mini-stirrer for mixing quantities up to 250 ml.

- Durable, brushless motor
- Continuously adjustable speed range
- High magnetic adhesion
- Portable operation with long operating time (8-12 h)
- Short charging time (2-3 h)
- Standard replaceable AA rechargeable batteries
- Optional power mode:
 - a) Mains-free with standard batteries
 - b) With supplied mains adapter (without batteries)
 - c) Combined mains/battery operation (with batteries fitted)

Accessories: TRIKA® 25 Magnetic stirring bar, RSE Stirring bar remover, IKAFLON® 15 Magnetic stirring bar, IKAFLON® 25 Magnetic stirring bar, IKAFLON® 10 Magnetic stirring bar, IKAFLON® 20 Magnetic stirring bar, IKAFLON® 30 Magnetic stirring bar

| Technical Data | |
|--|----------------|
| Number of stirring positions | 1 |
| Speed range [1/min] | 300 - 1800 |
| Set-up plate material | synthetic (PP) |
| Protection class according to DIN EN 60529 | IP 20 |
| Voltage [V] | 100 - 240 |
| Frequency [Hz] | 50/60 |
| DC Voltage [V=] | 9 |
| Power input [mA] | 1000 |
| Ident. No. | 3381300 |